



# MANUAL PENGGUNA

## USER MANUAL

Portable Mesh Nebulizer / Nebulizer Mesh Mudah Alih  
Model: LT-N400

Bahasa Malaysia + English

Prepared for Lionew Sdn. Bhd. | Dr Isla Brand

For home care use and medical institution use. Read carefully before use.



## KANDUNGAN / CONTENTS

N400

|    |                                                                   |
|----|-------------------------------------------------------------------|
| 1  | Keselamatan / Safety Guide                                        |
| 2  | Gambaran Produk / Product Overview                                |
| 3  | Spesifikasi & Kandungan Pakej / Specifications & Package Contents |
| 4  | Rupa Produk / Product Appearance                                  |
| 5  | Bekalan Kuasa / Power Supply                                      |
| 6  | Cara Penggunaan / How to Use                                      |
| 7  | Pembersihan & Disinfeksi / Cleaning & Disinfection                |
| 8  | Penyelesaian Masalah / Troubleshooting                            |
| 9  | EMC, Simbol, Pelupusan / EMC, Symbols, Disposal                   |
| 10 | Waranti & Maklumat Hubungan / Warranty & Contact                  |





## PANDUAN KESELAMATAN / SAFETY GUIDE

N400

### Bahasa Malaysia

Baca manual ini dengan teliti sebelum menggunakan produk.  
Produk ini adalah peranti perubatan aktif bukan implan untuk kegunaan nebulisasi.  
Pesakit atau pengguna awam perlu mendapatkan nasihat profesional kesihatan sebelum penggunaan.  
Gunakan ubat, dos dan cara penggunaan seperti diarahkan oleh doktor atau profesional kesihatan.  
Jangan gunakan ubat yang tidak dipreskripsikan, cecair pekat, cecair bersuspensi atau cecair berkepekatan tinggi.  
Jangan kongsi aksesori nebulizer dengan orang lain kerana risiko jangkitan silang.  
Kanak-kanak perlu menggunakan produk di bawah pengawasan orang dewasa.

### English

Read this manual carefully before using the product.  
This product is an active, non-implantable medical device for nebulization use.  
Patients or lay users should consult a healthcare professional before use.  
Use the medicine, dosage and operation method as advised by a doctor or healthcare professional.  
Do not use non-prescribed medicine, highly viscous liquid, suspended liquid or high-concentration liquid.  
Do not share nebulizer accessories with other users due to cross-infection risk.  
Children must use the product under adult supervision.

## AMARAN & BATASAN / WARNINGS & LIMITATIONS

N400

### Bahasa Malaysia

Jangan gunakan produk berhampiran sistem anestesia, ventilator, oksigen perubatan atau persekitaran MRI.  
Jangan jatuhkan, hentak kuat, ubah suai atau buka bahagian produk.  
Jangan rendam unit utama atau penyesuai kuasa di dalam air.  
Jangan sentuh mesh nebulizer dengan kapas, berus, pin atau objek tajam.  
Jangan bawa atau simpan peranti semasa masih mempunyai cecair ubat di dalam cawan ubat.  
Simpan produk jauh daripada cahaya matahari langsung, suhu melampau, haba, habuk dan lint.  
Gunakan hanya kabel/penyesuai kuasa yang mematuhi IEC 60601-1 dan sesuai dengan spesifikasi produk.

### English

Do not use the product near anaesthesia systems, ventilators, medical oxygen or MRI environments.  
Do not drop, strike, modify, disassemble or repair the product by yourself.  
Do not immerse the main unit or power adapter in water.  
Do not touch the nebulizer mesh with cotton swabs, brushes, pins or sharp objects.  
Do not carry or store the device while liquid medicine remains in the medication cup.  
Keep the product away from direct sunlight, extreme temperature, heat, dust and lint.  
Use only cables/adapters that comply with IEC 60601-1 and match the product specifications.

Important / Penting: Serious incidents related to the device should be reported to the manufacturer and the competent authority where the user or patient is established.

## GAMBARAN PRODUK / PRODUCT OVERVIEW



BM

**Prinsip kerja:** Nebulizer menggunakan getaran ultrasonik dan mesh mikro berpori untuk menukarkan cecair ubat kepada zarah aerosol halus.

**Tujuan penggunaan:**

Untuk terapi penyedutan larutan ubat yang dipreskripsikan kepada pesakit dewasa dan kanak-kanak berumur 2 tahun ke atas.

**Persekitaran:**

Sesuai untuk institusi perubatan dan penjagaan di rumah.

N400

EN

**Working principle:** The nebulizer uses ultrasonic vibration and a micro-porous mesh to convert liquid medicine into fine aerosol particles.

**Intended purpose:**

For inhalation therapy of physician-prescribed solutions for adult and paediatric patients aged 2 years and above.

**Environment:**

Suitable for medical institutions and home care use.

Contraindication / Kontraindikasi: Do not use if allergic to aerosolized medicines. Do not use Pentamidine powder/suspension. For asthma or special conditions, use only under healthcare professional advice.

## SPESIFIKASI PRODUK / PRODUCT SPECIFICATIONS

N400


| Specifications / Spesifikasi             |                                                         |
|------------------------------------------|---------------------------------------------------------|
| Perkara / Item                           | Maklumat / Information                                  |
| Model / Model                            | LT-N400                                                 |
| Nama produk / Product name               | Ultrasonic Mesh Nebulizer                               |
| Kaedah kerja / Working method            | Piezoelectric resonance                                 |
| Frekuensi resonan / Resonant frequency   | 110 kHz $\pm$ 10%                                       |
| Kadar nebulisasi / Nebulization rate     | $\geq$ 0.20 mL/min                                      |
| Saiz zarah median / Median particle size | 3.0 $\mu$ m $\pm$ 25%                                   |
| Taburan zarah / Particle distribution    | 1 $\mu$ m–5 $\mu$ m >50%                                |
| Kapasiti cawan ubat / Cup capacity       | Max. 8 mL; Min. 0.5 mL                                  |
| Bekalan kuasa / Power supply             | Built-in 3.7V lithium battery or DC 5V – 1A adapter     |
| Masa kerja / Working time                | 10 minutes/session; at least 30 minutes per full charge |
| Dimensi / Dimension                      | 48 x 43 x 107 mm                                        |
| Berat / Weight                           | Approx. 87 g                                            |

## KANDUNGAN PAKEJ & RUPA PRODUK / PACKAGE & APPEARANCE

N400

### Package contents / Kandungan pakej

| Item                            | Qty / Kuantiti |
|---------------------------------|----------------|
| Main unit / Unit utama          | 1              |
| Medication cup / Cawan ubat     | 1              |
| USB cable / Kabel USB           | 1              |
| Adult mask / Topeng dewasa      | 1*             |
| Child mask / Topeng kanak-kanak | 1*             |
| User manual / Manual pengguna   | 1              |




| Medication cup / Cawan ubat      | Medication cup button / Butang cawan ubat |
|----------------------------------|-------------------------------------------|
| Nebulizer sheet / Mesh nebulizer | USB interface / Port USB                  |
| Power button / Butang kuasa      | Warning indicator / Indikator amaran      |

## PENGEKASAN BATERI / BATTERY INSTALLATION

N400

### Using external power supply

1. Connect one end of the USB cable to the Type-C port of the main unit and the other end to the USB power adapter.
2. Plug the USB adapter into the power outlet.
3. While charging, the indicator light flashes blue, and when fully charged, it remains steadily blue.
4. When fully charged, the lithium battery can work for at least 30 minutes. The machine will automatically shut down every 10 minutes for safety protection purposes.
5. After approximately 300 charge-discharge cycles, the battery capacity may decrease to below 80%.



### ⚠ Note:

- Be sure to use the power adapter that meets the requirements of the IEC 60601-1 standard.
- The nebulizer should not be used during charging.

### Bahasa Malaysia

Sambungkan kabel USB Type-C pada port utama dan penyesuai kuasa USB yang sesuai.

Pasangkan penyesuai USB ke soket kuasa.

Semasa pengecasan, lampu indikator berkelip biru; selepas penuh, lampu kekal biru.

Apabila dicas penuh, bateri litium boleh digunakan sekurang-kurangnya 30 minit.

Peranti akan mati secara automatik selepas 10 minit untuk perlindungan keselamatan.

Nebulizer tidak boleh digunakan semasa pengecasan.

### English

Connect the USB Type-C cable to the main unit and a suitable USB power adapter.

Plug the USB adapter into the power outlet.

During charging, the indicator flashes blue; when fully charged, it remains steady blue.

When fully charged, the lithium battery can operate for at least 30 minutes.

The device shuts off automatically after 10 minutes for safety protection.

The nebulizer must not be used during charging.




## SEBELUM GUNA & PASANG TOPENG / BEFORE USE & MASK SETUP

N400

se refer to the user manual for instructions on how to use the mask.  
component needs to be cleaned and dried before installation.

ach the mask to the medication cup.



e : Please select the correct mask.

- a) The large mask is suitable for adults, and the small mask is suitable for patients.
- b) The mask can be reused 20 times. If necessary, please contact our company to purchase it through regular channels such as drug stores and hospitals.
- c) To remove, follow the steps in reverse order.
- d) Compatible device: Huizhou Kaiyi Technology Co., Ltd. Manufacturer type: NC180801.

the liquid into the cup and fasten the lid to prevent the liquid from coming out.



Open the medication cup cover

② Fill the medication into the cup

③ Cover and fasten the lid

### Bahasa Malaysia

Pastikan semua komponen telah dibersihkan dan dikeringkan sebelum pemasangan.  
Pilih topeng yang sesuai: topeng besar untuk dewasa dan topeng kecil untuk kanak-kanak.  
Pasangkan topeng pada cawan ubat dengan kemas.  
Topeng boleh digunakan sehingga 20 kali. Tukar topeng jika retak, rosak atau berubah bentuk.  
Jika perlu, beli aksesori yang serasi daripada syarikat atau saluran pembelian yang diluluskan.  
Untuk membuka semula, lakukan langkah pemasangan secara terbalik.

### English

Ensure all components are cleaned and dried before installation.  
Select the correct mask: large mask for adults and small mask for paediatric users.  
Attach the mask securely to the medication cup.  
The mask may be used up to 20 times. Replace the mask if it is cracked, damaged or deformed.  
If needed, purchase compatible accessories from the company or approved purchase channel.  
To remove the mask, follow the installation steps in reverse order.

## ISI CECAIR UBAT / FILLING MEDICATION CUP

N400

### Bahasa Malaysia

Buka penutup cawan ubat.  
Masukkan cecair ubat yang dipreskripsikan ke dalam cawan ubat.  
Pastikan isipadu tidak melebihi 8 mL dan tidak kurang daripada 0.5 mL.  
Tutup dan kunci penutup dengan kemas untuk mengelakkan kebocoran.  
Jangan guna cecair pekat, melekit, bersuspensi atau berkepekatan tinggi kerana boleh menyumbat mesh.  
Jika ubat tertumpah, lap dengan kain kasa bersih dengan segera.

### English

Open the medication cup cover.  
Pour the prescribed liquid medicine into the medication cup.  
Make sure the volume does not exceed 8 mL and is not less than 0.5 mL.  
Close and fasten the lid securely to prevent leakage.  
Do not use thick, sticky, suspended or high-concentration liquids because they may clog the mesh.  
If medicine is spilled, wipe it immediately with clean gauze.

### ⚠ Note : Please select the correct mask

- a) The large mask is suitable for adults,
- b) The mask can be reused 20 times. If purchase it through regular channels
- c) To remove, follow the steps in reverse
- d) Compatible device: Huizhou Kaiyi Te NC180801.

### 2. Pour the liquid into the cup and fasten the lid to prevent leaking out.



① Open the medication cup cover




② Fill the medication cup with the liquid

### ⚠ Note

## PROSES NEBULISASI / NEBULIZATION

N400

### 3. Nebulization



1. Press the power button for short duration, the indicator light (green) will be on.
2. After the automatic self-cleaning procedure, the green indicator remains steady and nebulization starts.
3. The nebulizer operates for 10 minutes per session and will shut off automatically. To continue using the device, please press the power button again to restart.
4. When the nebulization process is complete or if you wish to stop it manually, please press the power button to turn off the device.

### Bahasa Malaysia

Tekan butang kuasa sebentar untuk menghidupkan peranti; lampu indikator hijau menyala.  
Selepas pembersihan automatik, lampu hijau stabil dan nebulisasi bermula.  
Peranti beroperasi 10 minit bagi setiap sesi dan akan mati secara automatik.  
Tekan butang kuasa sekali lagi untuk meneruskan sesi atau menghentikan penggunaan secara manual.  
Pegang nebulizer dengan stabil supaya cecair ubat bersentuhan dengan mesh nebulizer.  
Jika prestasi kabus menurun, semak paras cecair dan kebersihan cawan ubat.

### English


Press the power button briefly to turn on the device; the green indicator lights up.  
After automatic cleaning, the green light remains steady and nebulization begins.  
The device operates for 10 minutes per session and shuts off automatically.  
Press the power button again to continue the session or stop manually.  
Hold the nebulizer steadily so the liquid medicine remains in contact with the nebulizer mesh.  
If the mist performance decreases, check the liquid level and cleanliness of the medication cup.

## PEMBERSIHAN & DISINFEKSI / CLEANING & DISINFECTION

N400

Do not share your nebulizer accessories, otherwise there is a risk of spreading infection and illness among different users. Please refer to the user manual for instructions on cleaning the mask.

### Cleaning the accessories upon first use and after each application.

- 
- ① Remove the medication cup from the main unit.
- 
- ② Pour out the remaining liquid.
- 
- ③ Clean the medication cup and mask with purified water not exceeding 40 °C.
- 
- ④ Dry the medication cup in a ventilated place, and then put it on the main unit.

### Cleaning the main unit

The main unit must be cleaned with a soft dry cloth and non-abrasive cleaners. Let it air dry on a clean paper towel.

#### Note:

1. Clean all parts of the nebulizer after each treatment to reduce the risk of infection.

### Bahasa Malaysia

Jangan kongsi aksesori nebulizer untuk mengurangkan risiko jangkitan silang.

Tanggalkan cawan ubat daripada unit utama selepas digunakan. Buang baki cecair ubat mengikut peraturan tempatan. Bersihkan cawan ubat dan mesh dengan air bersih yang tidak melebihi 40°C.

Keringkan cawan ubat di tempat berventilasi sebelum dipasang semula.

Bersihkan unit utama dengan kain lembut dan bahan pembersih tidak melelas; jangan rendam unit utama.

Bersihkan atau gantikan aksesori selepas setiap penggunaan.

### English

Do not share nebulizer accessories to reduce the risk of cross-infection.

Remove the medication cup from the main unit after use.

Dispose of the remaining liquid medicine according to local regulations.

Clean the medication cup and mesh with clean water not exceeding 40°C.

Dry the medication cup in a ventilated place before reinstalling it.

Clean the main unit with a soft cloth and non-abrasive cleaner; do not immerse the main unit.

Clean or replace accessories after each use.

## PENYELESAIAN MASALAH / TROUBLESHOOTING

N400

| Masalah / Failure                                                    | Punca dan tindakan / Possible cause and solution                                                                                                                                                                                                         |
|----------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Tidak berfungsi apabila dihidupkan / Does not work when powered on   | Cas semula bateri. Semak paras ubat dalam cawan. Pastikan cawan ubat dipasang dengan betul. / Recharge the battery. Check liquid level. Ensure medication cup is installed correctly.                                                                    |
| Kadar nebulisasi terlalu rendah / Nebulization rate is very low      | Pastikan cecair larut air dan tidak mengakis. Semak jika cecair mencukupi. Bersihkan atau ganti cawan ubat jika mesh tersumbat. / Ensure liquid is water-soluble and non-corrosive. Check liquid level. Clean or replace the cup if the mesh is clogged. |
| Indikator amaran kuning / Yellow warning indicator                   | Paras cecair ubat mungkin tidak mencukupi atau cawan ubat tidak dipasang dengan betul. / Liquid level may be insufficient or medication cup may not be installed correctly.                                                                              |
| Berhenti automatik semasa digunakan / Stops automatically during use | Cas bateri, isi semula ubat, pasang cawan ubat semula dan pegang peranti dengan stabil. / Recharge battery, refill medicine, reinstall the cup and hold steadily.                                                                                        |

Jika masalah masih berterusan, hentikan penggunaan dan hubungi pengedar / manufacturer. / If the problem persists, stop use and contact the distributor / manufacturer.

# SIMBOL, PELUPUSAN & EMC / SYMBOLS, DISPOSAL & EMC

N400

## Key symbols / Simbol utama

### 3. Explanation of Symbols

| S/N | Symbol                                                                              | Explanation                                                                                                                                     |
|-----|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| 1   |    | BF type application part                                                                                                                        |
| 2   |    | Note: refer to the attached file                                                                                                                |
| 3   |    | Refer to the manual                                                                                                                             |
| 4   |    | Serial number                                                                                                                                   |
| 5   | <b>IP22</b>                                                                         | Shell protection grade                                                                                                                          |
| 6   |    | Disposal of waste electrical and electronic equipment separately (follow local government regulations and recycling instructions for batteries) |
| 7   |    | Non-ionizing radiation                                                                                                                          |
| 8   |    | Date of manufacture                                                                                                                             |
| 9   |   | Manufacturer                                                                                                                                    |
| 10  |  | Use-by date                                                                                                                                     |
| 11  |  | Comply with the Regulation (EU) 2017/745                                                                                                        |
| 12  |  | Authorized representative in the European Union                                                                                                 |
| 13  |  | Medical device                                                                                                                                  |
| 14  |  | Unique device identifier                                                                                                                        |
| 15  |  | Single patient multiple use                                                                                                                     |

## EMC and disposal / EMC dan pelupusan

### BM

- Gunakan peranti mengikut maklumat EMC yang disediakan.
- Elakkan penggunaan berdekatan peralatan RF mudah alih, telefon bimbit, gelombang mikro dan sumber elektromagnet kuat.
- Jangan guna bersebelahan atau bertindan dengan peralatan lain kecuali telah disahkan berfungsi dengan normal.
- Lupuskan peranti, bateri dan aksesori mengikut peraturan tempatan.

### EN

- Use the device according to the EMC information provided.
- Avoid use near portable RF equipment, mobile phones, microwave ovens and strong electromagnetic sources.
- Do not use adjacent to or stacked with other equipment unless normal operation is verified.
- Dispose of the device, batteries and accessories according to local regulations.

Detailed EMC declaration tables are provided on the following pages. / Jadual deklarasi EMC terperinci disediakan pada halaman berikutnya.

## EMC DECLARATION / DEKLARASI EMC (1/6)

N400

Original Appendix 9.3 LT-N400 Manual — PDF page 17

BM: Gunakan peranti mengikut keperluan EMC. Elakkan penggunaan terlalu dekat dengan peralatan RF mudah alih atau susunan bertindan tanpa pengesahan operasi normal.

EN: Use the device according to EMC requirements. Avoid very close use to portable RF equipment or stacked use unless normal operation is verified.

### Chapter 8 / EMC

#### ⚠ Note:

- The instrument conforms to the requirements of IEC 60601- 1 - 2 standard for electromagnetic compatibility.
- The user shall install and use the EMC information provided in the random file.
- Portable and mobile RF communication equipment may affect the performance of the instrument. Avoid strong electromagnetic interference when using the device, such as being close to mobile phones, microwave ovens, etc.
- The guidance and manufacturer's declaration are detailed in the table below.

#### ⚠ Warnings:

- The instrument should not be close to or stacked with other equipment. If it must be close to or stacked, it should be observed and verified to operate normally under its configuration.
- In addition to the cables sold by the instrument manufacturers as spare parts for internal components, the use of other accessories and cables may result in increased emission or reduced immunity.

## EMC DECLARATION / DEKLARASI EMC (2/6)

N400

Original Appendix 9.3 LT-N400 Manual — PDF page 18

BM: Jadual ini menunjukkan deklarasi pelepasan elektromagnet dan imuniti awal bagi peranti LT-N400.

EN: This table shows the electromagnetic emissions and initial immunity declaration for the LT-N400 device.

### Chapter 8 / EMC

| Guidance and manufacture's declaration – electromagnetic emission                                                                                                                 |            |                                                                                                                                                                                                                                                                                                                |  |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| The product is intended for use in the electromagnetic environment specified below. The customer of the user of the product should assure that it is used in such an environment. |            |                                                                                                                                                                                                                                                                                                                |  |
| Emission Test                                                                                                                                                                     | Compliance | Electromagnetic environment – guidance                                                                                                                                                                                                                                                                         |  |
| Conducted and Radiated RF emissions CISPR 11                                                                                                                                      | Group 1    | The product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                                                                                                               |  |
| Conducted RF emissions CISPR 11                                                                                                                                                   | Class B    | The product is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes except for near active HF surgical equipment and the RF shielded room for magnetic resonance imaging. |  |
| Radiated RF emissions CISPR 11                                                                                                                                                    | Class B    |                                                                                                                                                                                                                                                                                                                |  |
| Harmonic emissions IEC 61000-3-2                                                                                                                                                  | Class A    |                                                                                                                                                                                                                                                                                                                |  |
| Voltage fluctuations and flicker emissions IEC 61000-3-3                                                                                                                          | Complies   |                                                                                                                                                                                                                                                                                                                |  |

| Guidance and manufacture's declaration – electromagnetic immunity                                                                                                                 |                                               |                                              |                                                                                                                                               |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|----------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment. |                                               |                                              |                                                                                                                                               |
| Immunity test                                                                                                                                                                     | IEC 60601-1-2 test level                      | Compliance level                             | Electromagnetic environment-guidance                                                                                                          |
| Electrostatic discharge IEC 61000-4-2                                                                                                                                             | ±8kV contact;<br>±2kV, ±4kV, ±8kV, ±15 kV air | ±8kV contact;<br>±2kV, ±4kV, ±8kV, ±15kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |



## EMC DECLARATION / DEKLARASI EMC (3/6)

N400

Original Appendix 9.3 LT-N400 Manual — PDF page 19

BM: Keperluan imuniti merangkumi transien cepat, lonjakan voltan, medan magnet frekuensi kuasa, penurunan voltan dan gangguan voltan.

EN: Immunity requirements cover fast transients, surges, power-frequency magnetic fields, voltage dips and voltage interruptions.

### Chapter 8 / EMC

|                                                                                                                                                                                                                              |                                                                                                                                     |                                                                                                                                  |                                                                                                                                                                                                                                                                                                     |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrical fast transients/bursts<br>IEC 61000-4-4                                                                                                                                                                           | ±2kV AC power supply lines;<br>±1kV DC power/Signal lines.<br>100 kHz repetition frequency                                          | ±2kV for AC power supply lines<br>100 kHz repetition frequency                                                                   | Mains power quality should be that of a typical commercial environment.                                                                                                                                                                                                                             |
| Surges<br>IEC 61000-4-5                                                                                                                                                                                                      | ±0.5kV, ±1kV lines to lines;<br>±0.5kV, ±1kV, ±2kV lines to earth                                                                   | ±0.5kV, ±1kV lines to lines                                                                                                      | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                 |
| Rated power frequency magnetic fields<br>IEC 61000-4-8                                                                                                                                                                       | 30A/m<br>50Hz or 60Hz                                                                                                               | 30A/m<br>50Hz or 60Hz                                                                                                            | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial environment.                                                                                                                                                                       |
| Voltage dips<br>IEC 61000-4-11                                                                                                                                                                                               | 0% UT, 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°;<br>0% UT, 1 cycle and 70% UT, 25/30 cycle<br>Single phase: at 0° | 0% UT, 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°; 0% UT, 1 cycle and 70% UT, 25/30 cycle<br>Single phase: at 0° | Mains power quality should be that of a typical commercial environment.<br>If the user of the product requires continued operation during power mains interruptions, it is recommended that the Electronic Stimulator for Pain Relief be powered from an uninterruptible power supply or a battery. |
| Voltage interruptions<br>IEC 61000-4-11                                                                                                                                                                                      | 0% UT, 250/300 cycle                                                                                                                | 0% UT, 250/300 cycle                                                                                                             |                                                                                                                                                                                                                                                                                                     |
| <p><b>⚠ Note:</b> UT is the a.c. mains voltage prior to application of the test level.<br/>E.g.: 25/30 means 25 periods at 50 Hz or 30 periods at 60 Hz.<br/>250/300 means 250 periods at 50 Hz or 300 periods at 60 Hz.</p> |                                                                                                                                     |                                                                                                                                  |                                                                                                                                                                                                                                                                                                     |

## EMC DECLARATION / DEKLARASI EMC (4/6)

N400

Original Appendix 9.3 LT-N400 Manual — PDF page 20

BM: Peralatan komunikasi RF mudah alih perlu dijauhkan daripada peranti dan kabel seperti dinyatakan dalam panduan EMC.

EN: Portable RF communications equipment should be kept away from the device and cables as stated in the EMC guidance.

### Chapter 8 / EMC

| Guidance and manufacture's declaration – electromagnetic immunity                                                                                                                                                                           |                                                                                                                                       |                                                                                                                                       |                                                                                                                                                                                                                |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.                                                           |                                                                                                                                       |                                                                                                                                       |                                                                                                                                                                                                                |
| Immunity Test                                                                                                                                                                                                                               | IEC 60601-1-2 test level                                                                                                              | Compliance Level                                                                                                                      | Electromagnetic Environment-Guidance                                                                                                                                                                           |
| Conducted disturbances induced by RF fields IEC 61000-4-6                                                                                                                                                                                   | 3Vrms in 0.15MHz – 80MHz; 6Vrms in ISM and amateur radio bands between 0.15MHz and 80MHz (Home healthcare environment) 80% AM at 1kHz | 3Vrms in 0.15MHz – 80MHz, 6Vrms in ISM and amateur radio bands between 0.15MHz and 80MHz (Home healthcare environment) 80% AM at 1kHz | Portable and mobile RF communications equipment should be used no closer to any part of these devices, including cables. Interference may occur in the vicinity of equipment marked with the following symbol. |
| Radiated RF EM fields IEC 61000-4-3                                                                                                                                                                                                         | 10V/m (Home healthcare environment), 80MHz – 2.7GHz 80% AM at 1kHz                                                                    | 10V/m (Home healthcare environment) 80MHz – 2.7GHz 80% AM at 1kHz                                                                     |                                                                                                                                                                                                                |
| Proximity magnetic fields IEC 61000-4-39                                                                                                                                                                                                    | 30 kHz, 8A/m<br>134,2 kHz, 65A/m<br>13,56 MHz, 7.5 A/m                                                                                | 30 kHz: 8A/m<br>134,2 kHz:65A/m<br>13,56 MHz: 7.5 A/m                                                                                 | /                                                                                                                                                                                                              |
| <b>Note 1:</b> At 80 MHz and 800 MHz, the higher frequency range applies.<br><b>Note 2:</b> These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people. |                                                                                                                                       |                                                                                                                                       |                                                                                                                                                                                                                |

# EMC DECLARATION / DEKLARASI EMC (5/6)

N400

Original Appendix 9.3 LT-N400 Manual — PDF page 21

BM: Nota teknis EMC berkaitan jalur ISM, kekuatan medan dan syarat ujian dikekalkan seperti manual asal.

EN: EMC technical notes on ISM bands, field strengths and test conditions are retained as in the original manual.

## Chapter 8 / EMC

- A The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.  
The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- B Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the These devices is used exceeds the applicable RF compliance level above, the These devices should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the These devices.
- C Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.
- D This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT.
- E The carrier shall be modulated using a 50 % duty cycle square wave signal.
- F r.m.s., before modulation is applied.

### Guidance and manufacture's declaration – electromagnetic immunity

| The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment. |                          |                  |                                      |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|------------------|--------------------------------------|
| Immunity Test                                                                                                                                                                     | IEC 60601-1-2 test level | Compliance Level | Electromagnetic Environment-Guidance |
| Proximity fields from RF wireless communications equipment<br>IEC 61000-4-3                                                                                                       | See the following table  | Complies         |                                      |

## EMC DECLARATION / DEKLARASI EMC (6/6)

N400

Original Appendix 9.3 LT-N400 Manual — PDF page 22

BM: Jadual medan proximity untuk peralatan komunikasi RF tanpa wayar dikekalkan sepenuhnya supaya tiada maklumat tercicir.

EN: The proximity fields table for RF wireless communications equipment is retained in full so no information is omitted.

Chapter 8

EMC

| Test frequency (MHz) | Band <sup>a)</sup> (MHz) | Service <sup>a)</sup>                                           | Modulation                                       | Immunity Test Level (V/m) |
|----------------------|--------------------------|-----------------------------------------------------------------|--------------------------------------------------|---------------------------|
| 385                  | 380 - 390                | TETRA 400                                                       | Pulse modulation <sup>b)</sup> 18 Hz             | 27                        |
| 450                  | 430 - 470                | GMRS 460, FRS 460                                               | FM <sup>c)</sup> ± 5 kHz deviation<br>1 kHz sine | 28                        |
| 710                  | 704 - 787                | LTE Band 13, 17                                                 | Pulse modulation <sup>b)</sup> 217 Hz            | 9                         |
| 745                  |                          |                                                                 |                                                  |                           |
| 780                  |                          |                                                                 |                                                  |                           |
| 810                  | 800 - 960                | GSM 800/900,<br>TETRA 800, iDEN 820,<br>CDMA 850,<br>LTE Band 5 | Pulse modulation <sup>b)</sup> 18 Hz             | 28                        |
| 870                  |                          |                                                                 |                                                  |                           |
| 930                  |                          |                                                                 |                                                  |                           |
| 1720                 | 1700 - 1990              | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS | Pulse modulation <sup>b)</sup> 217 Hz            | 28                        |
| 1845                 |                          |                                                                 |                                                  |                           |
| 1970                 |                          |                                                                 |                                                  |                           |
| 2450                 | 2400 - 2570              | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7            | Pulse modulation <sup>b)</sup> 217 Hz            | 28                        |
| 5240                 | 5100 - 5800              | WLAN 802.11 a/n                                                 | Pulse modulation <sup>b)</sup> 217 Hz            | 9                         |
| 5500                 |                          |                                                                 |                                                  |                           |
| 5785                 |                          |                                                                 |                                                  |                           |

**Note :** If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, the carrier may be pulse modulated using a 50 % duty cycle square wave signal at 18 Hz, while it does not represent actual modulation, it would be worst case.

19



## WARANTI & MAKLUMAT HUBUNGAN / WARRANTY & CONTACT

N400

### Bahasa Malaysia

Jika produk mengalami kerosakan, hubungi pengedar tempatan atau manufacturer.  
Unit utama mempunyai waranti 1 tahun dari tarikh jualan.  
Cawan ubat termasuk mesh nebulizer mempunyai waranti 3 bulan.  
Aksesori pakai habis seperti topeng tidak dilindungi waranti.  
Waranti tidak melindungi kerosakan akibat jatuh, hentakan, terkena air, penggunaan tidak betul, kemalangan atau pengubahsuaian struktur asal.  
Simpan rekod pembelian dan kad waranti untuk tujuan servis.

### English

If the product has a fault, contact the local dealer or manufacturer.  
The main unit has a 1-year warranty from the date of sale.  
The medication cup, including the nebulizer mesh, has a 3-month warranty.  
Consumable accessories such as masks are not covered by the warranty.  
Warranty does not cover damage caused by dropping, impact, water exposure, improper operation, accident or alteration of the original structure.  
Keep the purchase record and warranty card for service purposes.

### Contact Information / Maklumat Hubungan

LIONEW SDN. BHD. (202101040236 / 1440536-M)  
Unit B-12-13, Megan Avenue II, Jalan Yap Kwan Seng, 50450 Kuala Lumpur, Malaysia.  
Email: sales@aolon.my Tel: +603-30021329



## REKOD REVISI / REVISION RECORD

N400

| Document / Dokumen            | Manual Pengguna / User Manual – LT-N400 |
|-------------------------------|-----------------------------------------|
| Brand / Jenama                | Dr.Isla                                 |
| Company / Syarikat            | Lionew Sdn. Bhd.                        |
| Version / Versi               | 00                                      |
| Date / Tarikh                 | 12/8/2026                               |
| Prepared by / Disediakan oleh | Siti Fazlina Binti Abd Wahab            |
| Approved by / Diluluskan oleh | Sim Chi Song                            |

## KAD WARANTI / WARRANTY CARD

N400

### Bahasa Malaysia

Jika produk rosak, hubungi pengedar tempatan atau manufacturer.

Unit utama: waranti 1 tahun dari tarikh jualan.

Cawan ubat termasuk nebulizer sheet: waranti 3 bulan.

Caj pembaikan dikenakan selepas tempoh waranti tamat.

Aksesori pakai habis seperti topeng tidak dilindungi waranti.

Waranti tidak melindungi kerosakan akibat jatuh/hentakan, air/basah, operasi tidak betul, kemalangan atau ubah suai struktur asal.

### English

If the product has a fault, contact your local dealer or manufacturer.

Main unit: 1-year warranty from the date of sale.

Medicine cup including nebulizer sheet: 3-month warranty.

Repair fee may be charged after the warranty period expires.

Consumables such as masks are not covered by warranty.

Warranty excludes damage from dropping/impact, water/wet condition, improper operation, accident or alteration of original structure.

*Disclaimer: We are not responsible for any problems with the instrument if any parts are repaired by unqualified personnel.*

| Product         |                                                      |     | Model         |  |
|-----------------|------------------------------------------------------|-----|---------------|--|
| Production date |                                                      |     | Purchase Date |  |
| Invoice Number  |                                                      |     |               |  |
| Address         |                                                      |     |               |  |
| Contact         |                                                      | Tel |               |  |
| Service record  | Description                                          |     | Repairer      |  |
|                 |                                                      |     |               |  |
| Remarks         | Please present this card during machine maintenance. |     |               |  |